



At Zuverlasse, we go beyond supplying premium hygiene and disinfection products we redefine industry standards. Our products are meticulously engineered using cutting-edge technology, ensuring the highest levels of safety, superior performance, and seamless convenience for our customers but our commitment doesn't stop at delivering high-quality solutions. As a trusted leader in hygiene and disinfection, we are driven by a deeper purpose: to promote health, safety, and operational excellence across every environment we serve. With a relentless focus on innovation and expertise, Zuverlasse is setting new benchmarks in the industry empowering businesses with solutions they can truly rely on.

HYDRO PRO™ WIPES

HYDROGEN PEROXIDE DISINFECTANT WIPES

HAVING EN CLAIMS:

European Norms EN-17126: 2018 (E) Sporicidal for medical area .

European Norms EN-13727: 2012+A2:2015 (E) Bactericidal for medical area.

European Norms EN 1275:2005 (E) Fungicidal & Yeasticidal for medical area.

European Norms EN-14476:2013+A2:2019 (E) Virucidal for medical area.

European Norms EN-17387:2021 Bactericidal & Yeasticidal & Fungicidal for medical area.

European Norms EN-13624:2021 (E) Fungicidal & Yeasticidal for medical area.

European Norms EN-14563:2008 (E) Mycobactericidal & Tuberculocidal for medical area.

European Norms EN-14348: 2005 (E) Mycobactericidal for medical area.

European Norms EN-14561:2006 (E) Bactericidal for medical area.

European Norms EN-14562: 2006 (E) Fungicidal & yeasticidal for medical area.

European Norms EN-17111: 2018 (E) Virucidal for medical area.

European Norms EN-16777: 2018 (E) Virucidal for medical area, yeasticidal, virucidal for medical area.

European Norme EN-17848: 2023 Sporicidal .

European Norms EN-17846:2023 (E) Sporicidal for medical area.

European Norme EN-16615:2015 Bactericidal, Yeasticidal.

European Norms EN-21530:2004 (E) for material compatibility for most of the instrument surfaces in OT/ICU.

HYDRO PRO™ WIPES

HYDROGEN PEROXIDE DISINFECTANT WIPES



APPLICATION AREAS:

Patient-contact accessories, Gonioscopy lenses, Diagnostic contact lenses, Outer part of Flexible and Non-Flexible Video Laryngoscope Monitors and Holders, Intraoral scanners, Dental mirrors, Dental hand instrument handles, Dental curing light tips, Ultrasound probe holders, Probe cable connectors, Endoscope controls heads, Fiber-optics cables, Camera heads, Light guides, Endoscopic adapters, Air/Water syringes etc.

COMPOSITION:

Hydrogen Peroxide 6 %
Inactive Ingredients 94 %

USAGE:

Apply HYDRO PRO™ Wipes on the surface needed to be disinfected and wait till it gets dry Minimum 3 minutes required to achieve claims and 5 minutes for sporicidal claim.

COMPATABILITY:

Hydro Pro™ Wipes is compatible to the following material:

- Silicon
- Polyvinyl chloride (PVC)
- Polycarbonate (PC)
- Ceramic film
- Glass
- Polyoxymethylene (POM)
- Acrylonitrile butadiene styrene (ABS)
- Polypropylene (PP)
- Polymethyl methacrylate (PMMA)
- Stainless Steel 316
- Aluminium

Compatibility Method: The test item was evaluated for 14 days (±2 h) against multiple test specimens. Visual and tactile inspections confirmed that the following specimens showed no significant changes after exposure to Hydro Pro™, in accordance with the contact test described in ISO 21530:2004.



MDA
Technology



Bacteria/Virus



NON-
TOXIC



SKIN
FRIENDLY



DESCRIPTION:

HYDRO PRO™ Wipes disinfect and powerfully clean surfaces in one easy step with (Magnetic Dust Absorption) MDA Technology. Hydrogen peroxide technology eliminates 99.99% (percent) of bacteria and viruses. They are easy to kill pathogens on a broad range of surface, including vinyl, plastic, glass, stainless steel, sinks and more. The pre moistured, disposal wipes quickly remove soap scum, grease and grime throughout the workplace. HYDRO PRO™ Wipes are perfect for use on high-touch surfaces where germs can be easily transferred from person to person. HYDRO PRO™ Wipes ensure effective disinfection by staying wet for the recommended contact time.

- HYDRO PRO™ Wipes is a form of Escalated Hydrogen Peroxide (EHP™) Technology which is stabilized for prolonged efficacy.
- HYDRO PRO™ Wipes is for hospital, medical facilities and nursing homes. Primarily for use on nonporous medical surfaces, including baked enamel, chrome, glass, glazed ceramic, glazed porcelain, laminated surfaces, plastic, stainless steel and vinyl.
- HYDRO PRO™ Wipes provides a thorough clean with no harsh odours. Eliminates microorganisms and antibiotic-resistant organisms.
- HYDRO PRO™ Wipes are environment friendly wipes.
- HYDRO PRO™ Wipes are non - corrosive in nature with inbuilt rust inhibitors and corrosion inhibitors for enhanced protection.
- cleaning & disinfection at desired surface you can dispose wipes as per BMW laws.

PACK SIZE:

60 Wipes - SKU 7272
187 Wipes- SKU 7273
120 Wipes- SKU 7274

Zuverlāsse GmbH
Nürnberg, Germany.
management@zuverlassegmbh.de
contact@zuverlasse.com
sales@zuverlasse.com



www.zuverlasseindia.com



www.zuverlasse.solutions



Zuverlāsse Hygeine India Pvt. Ltd.
Plot No. F4, B/3, IDA, Kukatpally, Gandhi Nagar,
Hyderabad, Telangana – 500 037.

Toll Free Number: 1800-258-6457